



# Set up now:

## Be ready when you need care!

TytoCare gives you immediate access to a trusted provider from the comfort of home, saving you time and giving you peace of mind!

## Get started with your TytoHome in 3 simple steps:

- 1** **Download** the TytoCare app on your tablet or smartphone. **Click the Sign Up button** to create your account in the TytoCare app.
- 2** Follow the app instructions to pair your device using your Wi-Fi network. When prompted, enter the registration code: **Triad**
- 3** **Finish** the activation process by experiencing Tyto with Emma, your virtual guide. You'll learn how to use your TytoHome and complete your set-up.



**If you need help with set up or have questions, call 956-242-6075 to talk to a TytoCare support member, or call Frontier Member Services at 956-983-9272.**

If you haven't set up your TytoHome yet and need assistance,  
call 956-242-6075.

## Schedule a live visit:

1

**Schedule** a visit on the scheduling page here:

[Scheduling Page](#)

Make sure you select that you have a TytoCare device during the scheduling process.



2


**10 minutes** before the visit, get your kit ready and login to the TytoCare app.



3

Tap Start on the "Not Feeling Well?" card and select the "Start scheduled appointment" visit option.

You will be entered into the waiting room and your provider will start the visit at the scheduled time.



## Exam and forward

Do the exams **now** and send to your provider to review. No scheduled visit needed!

1

Login to the TytoCare app.

2

Tap Start on the "Not Feeling Well?" card and select "Send exams to your PCP" visit option.

3

Follow the prompts and complete the necessary exams.

4

Submit your exams to your provider for review. Your provider will be in touch within 1 business day.



# What's in the kit

TytoHome with  
**camera and  
thermometer**

Stethoscope  
adapter for **lung  
and heart exams**



Otoscope  
adapter for  
**ear exams**

Tongue depressor  
adapter for  
**throat exams**

